

# FlexiFrame

Stretch, Growth & Recovery Instrument



An accurate and flexible instrument for testing the stretch, growth and recovery properties of knitted and woven fabrics.

# FlexiFrame

## STRETCH, GROWTH & RECOVERY INSTRUMENT

We have created an instrument which is completely flexible while being accurate, standardised and calibrated at the same time. It is a comprehensive testing solution for **stretch, growth and recovery** where the use of a tensile tester would be prohibitive due to the time, usually a few hours.



### KEY BENEFITS

#### COMPLETE FLEXIBILITY

FlexiFrame allows testing of both woven and knitted fabrics on any station to a broad range of test methods.

#### STANDARDISED AND CALIBRATED

The Flexiframe can be calibrated to ISO 17025, giving you and your customers confidence in the repeatability and accuracy of your results.

#### OPTIONS TO SUIT EVERY LABORATORY

Supplied as either a portable 6 stations instrument, or wall mounted in multiples of 3, to suit your laboratory environment.

#### A RANGE OF ACCESSORIES

Accessories including weights and wire hangers provide the user with the tools they need to test to ASTM D3107, Arcadia AG31, Ralph Lauren and ASTM D2594.

#### INDEPENDENT TIMERS

Each station is equipped with a timer which can be numbered, allowing the user to take it away from the instrument to complete other tasks.

#### DIMENSIONAL STABILITY RULER

A ruler which measures stretch and shrinkage in both distance and percentage is an effective way to simplify the measurement of results.

#### 6 Station Wheel Base

Dimensions (mm)	Weight (kg)
Height: 1750	48.3*
Width: 1289	
Depth: 810	

\*Weight of instrument  
(without any attachments)

#### 3 Station Wall Mounted

Dimensions (mm)	Weight (kg)
Height: 1686	15.9*
Width: 553	
Depth: 205	

\*Weight of instrument  
(without any attachments)

### CONTACT US

#### For more information

web: [www.james-heal.co.uk/en/flexiframe](http://www.james-heal.co.uk/en/flexiframe)

tel: +44 (0) 1422 366 355

email: [sales@james-heal.co.uk](mailto:sales@james-heal.co.uk)